

# CitraCen

## Disinfectant Spray

**Tuberculocidal • Fungicidal • Bactericidal • Staphylocidal • Salmonellacidal • Pseudomonacidal • Virucidal\***

**CitraCen Disinfectant Spray** has been formulated for use on hard, nonporous, inanimate surfaces commonly found in facilities and areas such as: hospitals, nursing homes, medical and dental offices, blood banks, laboratories, emergency treatment facilities, veterinary clinics, animal shelters, kennels, ambulance and patient transfer vehicles, police and law enforcement vehicles, prisons and jails, prisoner holding areas, offices, funeral homes, restrooms, hotels, motels, restaurants, schools, day care facilities, hair care shops and salons.

**CitraCen Disinfectant Spray** is an effective disinfectant on hard, nonporous, inanimate surfaces and objects such as walls, floors, sinks, exterior bowl surfaces and exterior urinal surfaces, cabinets, tubs, showers, door handles, light switch covers, bed frames, wheelchairs, IV stands, over-bed tables, examination tables and chairs, bassinets, line carts, medicine carts, hamper, counter tops, telephones, waste containers, lockers, desks, tables, chairs, file cabinets, accessible surfaces of hydrotherapy bathing units, patient transfer lifts, clippers, combs, scissors, barber and salon chairs, trays, racks, cats, vehicle door handles, vehicle steering wheels, rubber floor mats, shovels, scrapers, nonporous rubber boots, nonporous rubber shoe covers, and nonporous rubber gloves.

Based on a contact time of 10 minutes at 25°C, this product is effective against *Staphylococcus aureus*, *Trichophyton mentagrophytes*, *Pseudomonas aeruginosa*, *Mycobacterium tuberculosis*, *Salmonella choleraesuis*, and the following \*viruses on pre-cleaned, hard, nonporous, inanimate surfaces: HIV-1 (AIDS virus), Herpes Simplex (Type 2), Influenza A2 Virus (Hong Kong), Avian influenza virus, Avian Bronchitis virus, Avian laryngotracheitis virus, Hog cholera virus, and Pseudorabies virus.

### DIRECTIONS FOR USE:

**As a disinfectant** – Thoroughly clean all surfaces with soap or detergent and rinse with water. Hold the can upright 6 to 8 inches from the surface. Saturate all surfaces with the recommended disinfecting solution, allowing to remain wet for 10 minutes before wiping, if wiping is required.

**As a deodorant:** Hold the can upright, then spray toward the center of the room for 3 to 5 seconds.

**As a mildew preventative:** Spray surface until moist.

### PACKAGE SIZES:

16 oz Aerosol Can (12 per case)

### Reorder No.

PHEN18539

### SPECIFICATIONS: for EPA Reg. No. 211-32

### EPA Est. No. 11694-KS-1

Product Designation:	A disinfectant and deodorant spray	
Form:	Aerosol	
Appearance:	Clear liquid	
Odor:	Citrus	
Active Ingredients:	o-phenylphenol	0.210%
	Ethyl alcohol	69.623%
		30.167%
Inert Ingredients:	Yes	
Biodegradable:	Yes	
Flash Point:	65°F (TCC)	
Freeze/Thaw Stability:	Will not freeze under normal ambient conditions	
Storage Stability:	Minimum shelf life, one year	
Corrosion Factor:	None, contains corrosion inhibitors	
Boiling Point:	165°F	
Product Base:	Alcohol	
Deodorizing Capacity:	Essential oils in combination with antimicrobial agents	
Foam Properties:	None	
Evaporation Rate:	Twice as fast as water	
Can Size:	20 ounces	
Net Contents:	16 ounces	



**CENTRAL SOLUTIONS, INC.**

401 Funston Road, P.O. Box 15276  
Kansas City, Kansas 66115  
Tel: (913) 621-6542 (800) 255-0262  
[www.centralsolutions.com](http://www.centralsolutions.com)

*Solutions for Personal and Environmental Hygiene*  
**Formulation • Development • Production**