

COOK EMERGENCY PNEUMOTHORAX SET



Used to treat simple and tension pneumothorax. Supplied sterile in peel-open packages. Intended for one-time use.




CATHETER INTRODUCER NEEDLE
Reinforced radiopaque FEP catheter and stainless steel, appropriate length needle



COOK CHEST DRAIN VALVE
Latex free




CONNECTING TUBE
Polyvinylchloride



ONE-WAY STOPCOCK
Plastic



SYRINGE



MOLNAR DISC WITH PULL TIE

SET ORDER NUMBER ¹	CATHETER		NEEDLE		CONNECTING TUBE	Global Product Number
	French Size	Length	Gage	Point	Length	
C-TPTS-8.5-6.0-FSNS	8.5	6.0 cm	15	lancet	30 cm	G13722

¹Set consists of items shown above and the following items: alcohol prep, povidone-iodine swabstick and porous transparent tape.