

THOMAS EMS TEMPERATURE CONTROL

▶▶ CLIMA-TECH™ CASE

**The Only PORTABLE TEMPERATURE CONTROLLED CASE
Designed Specifically for the EMS Market**



Everything comes securely packaged in a rugged Pelican case. Built with stainless steel hardware and padlock protectors for optimum security.



- Housed in a durable Pelican case
- Maintains 50 degrees off ambient temperature
- Automatically maintains + / - 1 degree from set temperature
- Retractable extension handle/polyurethane wheels
- Smart Sensor detects ambient temperature and predicts interior temperatures
- Will cool or heat up to 3 hours on battery power
- Automatically switches to battery power when removed from vehicle
- External digital temperature display
- Exterior Dimensions : 22" x 16.5" x 11"
- 12 Volt vehicle power or 110 volt with optional power supply

Inside the lid the contents are well organized. A double layer of elastic lined clear pockets assist with fast and secure placement of pharmaceuticals and other supplies.



The Clima-Tech's main storage compartment measures 14" W x 13" D x 6" H, providing adequate space for medications that require room temperature storage.



ITEM
#12348OR
#12348GR
#12348YL

DESCRIPTION
Clima-Tech Case, Orange
Clima-Tech Case, Gray
Clima-Tech Case, Yellow

PATENT PENDING

www.AllMed.net



888-633-6908