

LSP XTRA™ Backboard

Introducing the X-ray Translucent Rigid Affordable Backboard


The **XTRA™** backboard is one of the most significant developments in patient transport. For several years, many backboards have been designed on the basis of x-ray translucence (radio-transparency) instead of rigidity, at the expense of patient safety.

XTRA has been designed to provide the ultimate in x-ray translucence with the rigidity that is needed to protect the patient from further injury.

At 71" L x 15.75" W x 1.85" D, and weighing just 15 pounds, **XTRA** fits into standard storage areas and is easily carried to the patient.

The smooth surface allows for easy cleaning and disinfecting, while the construction of the **XTRA** prevents the intrusion of liquids into the board. **XTRA** will float an adult and has been tested to loads in excess of 1,000 lbs. without breaking.

Available with custom graphics in up to three colors, with no setup charges and no minimum order, your **XTRA** backboards can be customized with your logo and/or department or division.




Every body needs an XTRA Backboard!

Comparison of XTRA to Leading Competitors' Boards

Rigidity and Translucence

While many backboard manufacturers rightly claim either excellent x-ray translucence (which would place them high on the chart) or rigidity* (which would place them to the right on the chart), only XTRA has the perfect combination of x-ray translucence and rigidity.



*Based upon the deflection from horizontal of the center of the backboard with a 400 lb. distributed load, backboard supported at each end.

<h3 style="margin: 0;">SPECIFICATIONS</h3> <p>Length: 71" Width: 15.75" Depth: 1.85" Weight: 15.25 lbs. Capacity: 500 lbs. Material: High-density polymer</p>	<h3 style="margin: 0;">ORDERING INFORMATION</h3> <p>L812000 High Visibility Yellow, no pins L815000 Red, no pins L816000 Blue, no pins L818000 Orange, no pins L812014 High Visibility Yellow, 14 pins L815014 Red, 14 pins L816014 Blue, 14 pins L818014 Orange, 14 pins</p>
<p>Can be imprinted with custom logos in up to three colors, with no setup charges and no minimum order. Contact our Customer Service Department for assistance with ordering custom graphics.</p>	



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, Missouri 63110
 314/771-2400 • 800/444-3954 • Fax: 314/771-0650 • www.alliedhpi.com
 International offices in Europe, Asia, and South America

U.S.	Canada & Mexico	International
Toll Free: 800/225-4577	Toll Free: 800/446-0552	Telephone: 314/268-1683
Fax: 800/477-7701	Fax: 800/246-6201	Fax: 314/771-5183

©Allied Healthcare Products, Inc. Form No. 90-00-0004 8/04