

BUILT FOR LIFE™



20 Veterans Drive, Chicopee, MA 01022 | Tel.800.669.6900 | Fax.800.757.8585 | info@ironduck.com | www.ironduck.com

WEIGHT CAPACITY:

500 lbs. Weight Capacity

DIMENSIONS:

Standing: 72" x 16" x 2"

Folded: 36" x 16" x 4"

Weight: 19 lbs

DESIGN:

- 2 Rotomolded one piece seamless HDPE shell over non-toxic Eco-Friendly polyurethane foam*(1)
- 12 Over Sized Hand Holds
 - 12 Speed Pin slots
- 12 Pediatric strapping slots

WARRANTY:

Lifetime Warranty

RELATED PRODUCTS:

- Straps & Strapping Systems
- Sof Loc Head Immobilizer
- Backboard Sleeve
- Ultra Spac-Sav Carry Case

OPTIONS:

- P.P.E. Anti-Bac Additive
- Molded-In Graphics or
- Custom Screen Printing*(2)*(3)

STANDARD COLORS:

- Yellow, Blue, Orange, Red
- Lime Green, Olive Drab, Black



Ultra Spac-Sav

Product #: 35940
With Pins: 35940-P

Iron Duck's innovative folding plastic backboard can be backpacked for the long search and rescue haul or the urban stair climb and it stores in half the space, making it ideal for custom patient transport vehicles and quads. Featuring minimal deflection at weights in excess of 500 lbs, the Ultra Space-Sav is constructed with our exclusive stainless steel mechanical pivot and internal aluminum runners. The only rotomolded, folding backboard on the market, the durable HDPE shell does not reflect extreme temperatures ensuring superior patient comfort over metal folding stretchers. The subtle surface texture prevents slippage but is easy to clean and disinfect. Twelve oversized hand-holds each measuring 2" x 6" will accommodate gloved hands. 12 optional speed clips can be located in slots separate from hand-holds to facilitate accurate strapping and keep pins out of the way when lifting patients. Ships with 10 two inch reflective strips for added safety in low light situations – simply remove the backing from each strip and press into corresponding molded-in depressions around the edge of the backboard. Anti-bacterial agent can be molded in at additional cost. Lifetime Guarantee."

Now Available



- + P - Protects Patients
- + P - Protects EMS Professionals
- + E - Extends the Life of the Board

Anti-Staph, Anti-E Coli, Anti-Pneumonia... Anti-Bac



*(1) Kyoto & Montreal Compliant, *(2)Free on the Purchase of 25 or more Boards *(3)No Screen Printing on Black Backboards.
© Copyright 2009. Iron Duck, Built for Life, and Iron Duck logo are registered trademarks of Fleming Industries Inc. All Rights Reserved.