

The Ambu® Perfit ACE™ Collar



The Perfit ACE actually adjusts to all adult neck sizes with 16 settings (as well as the 4 standard adult size settings). The Perfit ACE incorporates latches, which allow the rescuer to size the collar according to the victim's exact neck size.



Main Features Include:

New Nasal Cannula Holders Molded into Collars

- Adjusts to adult size necks - ***16 settings***
- Stores completely flat
- Minimizes inventory
- Avoids using the wrong size collar
- Incorporates the Perfit numbering and color sizing system
- Horizontal sizing line
- Flip chin piece for easy airway access
- Radiolucent, CT, and MRI compatible
- Safety buttons allow for maximum reliability
- Nasal cannula holders located on collar body
- Part Number #281-000 (Perfit ACE - 30/case)

Ambu 

Ambu Inc. 6740 Baymeadow Drive, Glen Burnie, MD 21060
Phone: (800) 262-8462 Fax: (800) 262-8673
www.ambuUSA.com

The Ambu® Mini Perfit ACE® Collar

The Mini Perfit ACE actually adjusts to infant, pediatric, and small adult neck sizes with 12 settings as well as 3 standard sizes (infant, pediatric, & neckless). The Mini Perfit ACE incorporates latches, which allow the rescuer to size the collar according to the victim's exact neck size.



12 Adjustable Settings

Main Features Include:

The Mini Perfit ACE is perfect for SMALL ADULTS with thin neck sizes!

- Adjusts to infant, pediatric and small adult size necks – **12 settings**
- Stores completely flat
- Minimizes inventory
- Avoids using the wrong size collar
- Incorporates the Perfit numbering and color sizing system
- Horizontal sizing line
- Flip chin piece for easy airway access
- Radiolucent, CT, and MRI compatible
- Safety buttons allow for maximum reliability
- Nasal cannula holders located on collar body
- Part Number #281-106 (Mini Perfit ACE - 30/case)

Ambu

Ambu Inc. 6740 Baymeadow Drive, Glen Burnie, MD 21060
Phone: (800) 262-8462 Fax: (800) 262-8673 www.ambuUSA.com