

Introducing the NEW STANDARD in Oropharyngeal Evacuation...

THE SSCOR S3™ OROPHARYNGEAL EVACUATION TOOL

Scoops Chunks
Suctions Viscous Fluids
Sumps Pooling Fluids



**Replaces
The Finger Sweep.**



The SSCOR S3 is a rigid yet flexible large bore suction tip with a concave "scooping" surface which allows for removal of solids and other material that will not travel through a large bore suction tube.

**Assists in
Difficult Intubations.**



By placing the slider over the thumb port, the S3 can be briefly placed in the deep oropharynx to suction, freeing both hands to manage intubation.

Order #: 14369 -- each



888-633-6908 • www.AllMed.net • fax: 800-4ALLMED



LIT-44305 Rev. B



Emergo Europe
Molenstraat 15
2513 BH, The Hague
The Netherlands
Tel: (+31) 70 345 8570
Fax: (+31) 70 346 7299