

# CELOX™ GAUZE

## Rolled gauze coated with CELOX hemostatic agent

Introducing CELOX GAUZE, the latest in our innovative line of CELOX products. CELOX GAUZE combines our ground breaking CELOX agent with a fabric dressing to control life threatening blood loss directly at the source of bleeding.

CELOX GAUZE is the only fabric dressing available coated with the CELOX hemostatic agent. CELOX controls bleeding by forming a robust plug when red blood cells react with the hemostatic agent. The flexible gauze allows for easy and precise placement of the wound treatment overcoming the challenges typically associated with granular agents that must be poured into the wound.

CELOX's properties have been confirmed in wound models, where they have repeatedly demonstrated the ability to rapidly control major arterial bleeding (without cautery) from traumatic wounds.



### FEATURES AND BENEFITS

- Coated gauze allows for greater placement at source of bleeding
- Not affected by windy conditions
- Helps save lives by controlling bleeding within minutes
- Works in hypothermic conditions
- Works on anticoagulated blood
- No heat – No burning
- Sterile, waterproof packaging
- No special requirements
- Easy to remove and irrigate from wound

### TECHNICAL DATA

- One Gauze Strip: 3in x 10ft (7.6cm x 3.0M)

### PUBLICATIONS

**An Alternative Hemostatic Dressing: Comparison of CELOX, HemCon, and QuikClot;** Buddy G. Kozen, MD, LCDR, MC, USN, Sara J. Kircher, BS, RLAT, Jose Henao, MD, LCDR, MC, USN, Fermin S. Godinez, DO, Andrew S. Johnson, MD, CDR, MC, USN

### ORDERING INFO

Part Number	Description	Case Size
C12857	CELOX GAUZE Roll	25 pouches/case

ISO 13485 and 9001 Certified



**SAM MEDICAL PRODUCTS®**  
sammedical.com

customerservice@sammedical.com  
P.O. Box 3270 Tualatin, OR 97062  
800.818.4726 (USA) Tele  
503.639.5425 (USA) Fax

**SAM** MEDICAL  
PRODUCTS®

YOUR AUTHORIZED SAM MEDICAL PRODUCTS DEALER



Please call us at 888-633-6908 or visit our website, [www.AllMed.net](http://www.AllMed.net).