

C-A-T®

COMBAT APPLICATION TOURNIQUET®

PN: 30-0001 (Black) NSN: 6515-01-537-2611 (Black) PN: 30-0023 (Orange) PN: 30-0033 (Trainer Blue)

*"The Best
Prehospital Tourniquet"...*

February 2008 Supplement of The Journal of Trauma

*Over 4,000,000 sold to the
Department of Defense alone...*

The C-A-T® is a true one-handed tourniquet proven to be 100% effective by the U.S. Army's Institute of Surgical Research. Tests proved that the C-A-T® completely occluded blood flow of an extremity in the event of a traumatic wound with significant hemorrhage.

The C-A-T® utilizes a durable windlass system with a free-moving internal band providing true circumferential pressure to the extremity. Once adequately tightened, bleeding will cease and the windlass is locked into place. A Velcro® strap is then applied, securing the windlass to maintain pressure during casualty evacuation. The C-A-T®'s unique design avoids the use of screws and clips which can become difficult to operate under combat/survival stress or where fine motor skills are compromised.

- Awarded One of the "Top 10 Greatest Inventions" by the U.S. Army
- Small and Cost-Effective
- Designed to Perform in All Weather Conditions
- Proven to be 100% Effective in Occluding Blood Flow in both Upper & Lower Extremities by the U.S. Army's Institute of Surgical Research
- True, One-Handed Application



Laboratory Tested - Combat Proven®..

NORTH AMERICAN RESCUE™
INCORPORATED

www.NARescue.com 888.689.6277